




VYGON

Value Life

Enteral Nutrition
Nutrisafe²® Product Catalogue



The proven, safe, small and accurate enteral feeding system for neonates

Nutrisafe2®

The safe, small and most accurate nasogastric feeding product range for neonates

We have worked with leading clinicians in the field of enteral feeding to develop a complete range of nasogastric feeding tubes and accessories, allowing you to deliver the safest care for your patients.

- Small, lightweight, smooth connector surfaces specifically designed with neonates in mind
- Nutrisafe2 syringes are intended for oral and enteral application
- Delivers low dose accuracy of medication, without changing current policy and procedures
- Unique, proven design with 11 years global experience and 50 million connections globally every year
- Zero reported misconnections with any other small bore connector designs
- Accessories facilitate the preparation, transfer and administration of enteral fluids, providing a complete enteral system for neonates.

Contents


Enteral Feeding Tubes	3-5
Oral/enteral syringes	6
Accessories	7-11



Enteral Feeding Tubes

Nutrisafe2 feeding tubes are available in PVC, PUR and silicone for enteral nutrition of neonates and children.

- Features excellent radiopacity and a closed distal tip with two lateral eyelets
- The safe, small lightweight Nutrisafe2 connector reduces risk of rhinitis and oesophageal erosion compared to larger and heavier connectors
- Renowned soft textured feeding tubes
- Easy to insert
- 40cm, 50cm and 75cm length tubes include markings every centimetre from 5cm to 35cm. The 125cm length tubes include markings every centimetre from 5cm to 70cm for post-pyloric feeding
- Colour-coded tethered caps indicate the size of the feeding tube.



Enteral Feeding Tubes

Ordering Information

DEHP-free PVC Feeding Tube - For short-term use


Length (cm)	Size Fr (Outer Diameter)						Unit of Sale
	4Fr	5Fr	6Fr	8Fr	10Fr	12Fr	
40	0364042	0364052	0364062	0364082	N/A	N/A	25
50	0361042	0361052	0361062	0361082	0361102	N/A	50
75	0363042	0363052	0363062	0363082	0363102	0363122	25
125	0362042	0362052	0362062	0362082	0362102	0362122	25



0361042

Polyurethane Feeding Tube - For use up to four weeks

Length (cm)	Size Fr (Outer Diameter)						Unit of Sale
	4Fr	5Fr	6Fr	8Fr	10Fr	12Fr	
50	1361042	1361052	1361062	1361082	1361102	N/A	50
75	1363042	1363052	1363062	1363082	1363102	1363122	25
125	1362042	1362052	1362062	1362082	1362102	1362122	25



1361042

Silicone Feeding Tube - For use up to four weeks

Length (cm)	Size Fr (Outer Diameter)					Unit of Sale
	4Fr	5Fr	6Fr	8Fr	10Fr	
50	2331042	2331052	N/A	N/A	N/A	10
125	2332042	N/A	2332062	2332082	2332102	10
125 (W)	N/A	N/A	2395062	2395082	N/A	10



2331042

Silicone Duodenal weighted tip Feeding tube - For use up to four weeks

Length (cm)	Size Fr (Outer Diameter)		Unit of Sale
	6Fr	8Fr	
125	2395062	2395082	10

Ordering Information

DEHP-free PVC Salem Type Dual Flow Gastric Tubes

PVC dual-flow gastric tube for continuous gastric aspiration with Nutrisafe2 connections. 60cm with 5cm markings from 5cm to 45cm.

Length (cm)	Size Fr (Outer Diameter)		Unit of Sale
	6Fr	8Fr	
60	0340.062	0340.082	10



Oral and Enteral Syringes

Single-use and reusable Nutrisafe2 syringes

Ordering Information

Vygon Code	Product Description	Volume (ml)	Graduations (ml)	Connection	Unit of Sale																																			
I015013	Single-use Syringe	1	0.01	Centred	160																																			
I015023	Single-use Syringe	2.5	0.01	Centred	130																																			
I015053	Single-use Syringe	5	0.01	Centred	100																																			
I015103	Single-use Syringe	10	1	Offset	100																																			
I015213	Single-use Syringe	20	1	Offset	100																																			
I015602	Single-use Syringe	60	1	Offset </tr <tr> <td>I012.01</td> <td>Reusable Syringe</td> <td>1</td> <td>0.01</td> <td>Centred</td> <td>100</td> </tr> <tr> <td>I012.02</td> <td>Reusable Syringe</td> <td>2.5</td> <td>0.1</td> <td>Centred</td> <td>100</td> </tr> <tr> <td>I012.05</td> <td>Reusable Syringe</td> <td>5</td> <td>0.2</td> <td>Centred</td> <td>100</td> </tr> <tr> <td>I012.10</td> <td>Reusable Syringe</td> <td>10</td> <td>1</td> <td>Centred</td> <td>100</td> </tr> <tr> <td>I012.20</td> <td>Reusable Syringe</td> <td>20</td> <td>1</td> <td>Centred</td> <td>50</td> </tr> <tr> <td>I012.60</td> <td>Reusable Syringe</td> <td>60</td> <td>1</td> <td>Offset</td> <td>25</td> </tr>	I012.01	Reusable Syringe	1	0.01	Centred	100	I012.02	Reusable Syringe	2.5	0.1	Centred	100	I012.05	Reusable Syringe	5	0.2	Centred	100	I012.10	Reusable Syringe	10	1	Centred	100	I012.20	Reusable Syringe	20	1	Centred	50	I012.60	Reusable Syringe	60	1	Offset	25
I012.01	Reusable Syringe	1	0.01	Centred	100																																			
I012.02	Reusable Syringe	2.5	0.1	Centred	100																																			
I012.05	Reusable Syringe	5	0.2	Centred	100																																			
I012.10	Reusable Syringe	10	1	Centred	100																																			
I012.20	Reusable Syringe	20	1	Centred	50																																			
I012.60	Reusable Syringe	60	1	Offset	25																																			

Offset syringe tip on large syringe sizes



Positioning the syringe in a syringe driver with the tip above the centre point ensures nutrient rich fats that are administered first upon commencement of feed.

Smooth tip, smaller sizes

Small, smooth, atraumatic surfaces for oral application.



I015 range

Accessories

Drawing-up straw - 817154, 817052

The Nutrisafe2 sampling straws are connected to a Nutrisafe2 syringe for drug sampling (code 817.052) or to fill the syringe with milk/milk formula (other codes).

Some sampling devices have a smaller bore for non-viscous solutions. Other sampling devices have a larger bore for fortified milk.



Syringe cap single and 8 pack - 828.01, 828.08

Purple coloured obturator for Nutrisafe2 syringes, preserving enteral fluids prior to administration.



Extension Tubing - 0368.152, 0368.1521

PVC extension tubing provides extended connections to pumps, allowing greater flexibility when used at the cot side.



Vented caps - 818.22, 818.62

Vented cap to be fitted on Nutrisafe2 syringes (20ml and 60ml) to enable gravity feeding and protecting enteral feed from contamination.



Sampling needle - 817.002

Designed to aspirate drugs from an ampoule, which must be administered orally or via a Nutrisafe2 feeding tube.



Filter Straw - 817.105

Filter straw with Nutrisafe2 connection and filter (5µm). Designed to aspirate fluids from glass vials and ampoules whilst removing particulate matter.



Feeding tube obturator - 829.01

Purple-coloured obturator cap for Nutrisafe2 feeding tubes.



Spike - 0281062

Spike featuring (at the proximal end) a male Nutrisafe2 connection with cap. When connected to a Nutrisafe2 syringe, this device allows drug samplings from fluid reservoirs.



Drip stand holder - 367011

Holder for 20ml and 60ml Nutrisafe2 syringes for gravity feeding.



Three-way stopcock - 580201

Allows secondary connection when intermittent medication must be given via a Nutrisafe2 feeding tube.



Vial access caps - 821.012,022

Available for vial necks of 14mm or 20mm diameter for the purposes of drug sampling without the use of needles.



Male to male connector - 368.22

Male to male Nutrisafe2 connector to transfer drugs contained in a Nutrisafe2 syringe into another Nutrisafe2 syringe, or to transfer vitamins contained in a Nutrisafe2 syringe into the syringe containing the baby's milk.



Nutrisafe2/ENFit adaptor - 0368.72

The Nutrisafe2/ENFit adaptor enables the connection between an administration set equipped with a female ENFit connector and a Nutrisafe2 feeding tube.



Nutrisafe2/ENFit adaptor - 0VEACC510

The Nutrisafe2/ENFit adaptor allows the connection of Nutrisafe2 sets for pumps and Nutrisafe2 syringes to ENFit feeding tubes.



Neonatal Grip-Lok™ - 58041 I

Fixation device with hydrocolloid adhesive for enteral feeding tubes and nasal cannulae.

The device features the Grip-Lok technology that allows it to re-open it in order to gain access to the tube to clean or reposition it.

The Grip-Lok fixation device can be cut in half in order to reduce its size and adapt it to the patient.



Taper connector - 368.12

This taper connector with Nutrisafe2 male connection is used when it is necessary to administer medications with a Nutrisafe2 syringe into a gastrostomy device.



Funnel adaptor - 368.02

To connect onto funnel tubing.



Y-connector - 368.42

The Y-connector allows secondary connection when medication must be given with feed.



Battle cap for sampling - 820.20, 820.22 and 820.28

For Sampling for 20mm, 22mm and 28mm vials.



Gastric aspirator - 534.082

A collection trap volume 25ml, with a cap that ensures a tight seal of the specimen sampling.

Also features a tube connecting the device to the suction source and a funnel connector.

The tube has a female Nutrisafe2 hub to connect the device to the Nutrisafe2 feeding tube.



Ordering Information

Vygon Code	Product Description	Length (cm)	Bore Size	Unit of Sale
0817052	Drawing-up Straw	5	Standard	50
0817154	Drawing-up Straw	15	Fine	50
0817152	Drawing-up Straw	15	Standard	100
0817153	Drawing-up Straw	15	Large	50
0817202	Drawing-up Straw	20	Standard	50
0817203	Drawing-up Straw	20	Large	50
828.01	Syringe Cap Individually Packed			100
828.08	Syringe Cap in Packs of 8			100
0368032	Extension Tube	30	Standard	50
0368102	Extension Tube	100	Standard	50
0368152	Extension Tube	150	Standard	50
0369152	Extension Tube	150	Large	25
368.1021	Extension Tube with Clamp	100	Standard	50
368.1521	Extension Tube with Clamp	150	Standard	50
818.22	Vented Cap for 20ml Syringe			50
818.62	Vented Cap for 60ml Syringe			50
817.002	Sampling Needle			100
817.105	Sampling Device with Filter			50
829.01	Feeding Tube Obturator			100
0281062	Spike with Nutrisafe2			10
367011	Drip Stand Syringe Holder 20ml & 60ml Syringes			5
5802.01	Three-way Stopcock			50
821.012	Vial Access Cap with Nutrisafe2 Connection 14mm Male			50
821.022	Vial Access Cap with Nutrisafe2 Connection 20mm Male			50
0368.22	Male to male transfer device			30
0368.72	Nutrisafe2/ENFit adaptor			100
0VEACC510	Nutrisafe2/ENFit adaptor			300
580411	Neonatal Grip-Lok™ Fixation			50
368.12	Taper Connection			30
368.02	Funnel Adaptor			30
368.42	Y-connector			50
820.20	Bottle Cap for Sampling for 20mm Vials			20
820.22	Bottle Cap for Sampling for 22mm Vials			20
820.28	Bottle Cap for Sampling for 28mm Vials			20
534.082	Gastric Aspirator			10

For further information, please contact: vygon@vygon.co.uk

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon (UK) Ltd • The Pierre Simonet Building • V Park • Gateway North
Latham Road • Swindon • Wiltshire • SN25 4DL
Tel: 01793 748800 • Fax: 01793 748899 • Twitter: @vygonuk
Web: www.vygon.co.uk Content correct 03/2017 Code: DXJB0062002 v3
Copyright Vygon (UK) Ltd 2017



www.vygon.co.uk